

BRAINPET - A NEW INSERT FOR PET-MR PROTOTYPE

21.08.2018 | DAVID ARUTINOV



Member of the Helmholtz Association





CONTENT



- INTRODUCTION
- PET REQUIREMENTS FOR MR
- CONCEPT OF THE PET INSERT
- SUMMARY & STATUS






INTRODUCTION



MAIN IDEA

- The idea of the project is to mix two modalities, PET and MR inside a single device.
- Synergy of the two modalities brings simultaneous, non-invasive diagnostics to the next level as it combines morphological and functional imaging
- A completely new, high-resolution, and highly sensitive PET insert must be designed to operate in combination with UHF MR with the field up to 9.4T







INTRODUCTION



Positron

γ photon (511 keV)

WORKING PRINCIPLES OF TWO MODALITIES



- A body acts as a radio transmitter in a very high magnetic field if it is exposed to an EM wave.
- An image is obtained from the "echo radio transmission" of your body.
- Detection: mostly regions including fat and water
- A trace amounts of short-lived radioactive labeled molecules are injected into bloodstream

PET

γ photon (511 keV)

Electron

Radioactive

Radioactive decay can reveals dynamics of metabolic processes

Problems for electronics:

- Very high static magnetic fields (up to 9.4T)
- Additional dynamic magnetic field
- Very low EM signal



Problems for electronics:

- Complex / fast readout electronics
- Large amount of data to be processed



PET REQUIREMENTS FOR MR



REQUIREMENTS OF THE PET INSERT

Restrictions on the PET electronics:

- Works inside static Bo magnetic fields up to a 9.4T.

- ... fast changing gradient field of several tens of militesla.

- ... accompanying RF environment

Avoid using ferrites, mechanical cooling FANs, vacuum tubes, etc.

No real countermeasures are available, avoid large areas of conductive materials

restriction on a choice of

Proper shielding, but should be transparent to gradient field

Applies certain


data transmission

Requirements for PET Insert from the MR system:

should not emit certain frequencies and harmonics thereof.

 Parts close to the MR bore should contain as low amount of magnetic material as possible (at best, after 5 Gauss line).

- PET insert should be disinfectable, human safe (no current leakage), total weight is limited, etc.




safe (cm)

21 August 2018




DETECTORS



- PET module comprised of detector stacks
- Communication via optical link
- Air cooled
- Carbon shielded [2]






The detector stack based on dSiPM [2].

- 1. Scintillator
- 2. Philips dSiPM sensor module
- 3. FPGA based readout module
- 4. Powering





DETECTOR RING AND THE GANTRY










DATA READOUT AND CLOCK DISTRIBUTION



PCIe Gen.3 **Host Inter**face Card

PCle Gen.3

Server










DATA PROCESSING: MODE I



All the data is transmitted over the backend Interface and stored into an external memory (a. k. a. **debugging mode**)


Total memory bandwidth: ~ 6GB/s







DATA PROCESSING: MODE II




Different filters and algorithms reduce the data, and only part required for further analysis and image reconstruction gets transmitted over the backend Interface and stored in the external memory (server).







POWER SUPPLY



SUMMARY AND STATUS



- Development of the electronics for high magnetic fields requires a lot of care, since conventional designs might not be applicable. The literature is scarse, not much experienced people either.
- Our current concept went through several redesigns and we will have a final version soon.
- In order to eliminate other risks and better focus on project complexity we will try to base most of the hardware on already existing products on the market.

Thanks You!





REFERENCES



- [1] D. W. Harberts and M. van Helvoort, "MRI image distortion due to magnetic materials in medical i mplants," 2015 IEEE International Symposium on Electromagnetic Compatibility (EMC), Dresden, 2015, pp. 1463-1466.
- [2] Weissler, Björn, et al. "Design concept of world's first preclinical PET/MR insert with fully digital silicon photomultiplier technology." Nuclear Science Symposium and Medical Imaging Conference (NSS/MIC), 2012 IEEE. IEEE, 2012.



